- Each file has 3 composite spectrum plots for each precursor loss composite type, which are all intensity vs. m/z loss (found from subtracting measured m/z from precursor m/z)
 - 1. Composite spectrum zoomed in around water loss $z=2 (\pm 30 \text{ bins around theoretical m/z})$
 - 2. Composite spectrum zoomed in around phosphoric acid loss z=2 (± 30 bins around theoretical m/z)
 - 3. Full neutral loss composite spectrum, which measures from m/z of 0 to 100
- Plot details
 - o Blue curve is the composite spectrum from the file
 - o Red curve is the Gaussian fitted peak only present if a fit is found for the theoretical m/z
 - Vertical line (orange/green) is at the theoretical m/z of water loss z=2 (plot (1)) and phosphoric acid z=2 (plot (2)) orange if the peak is not significant, turns green if the fit is significant (criteria for water z=2 is >50 n_spectra, for phosphoric_acid z=2 is >50 n_spectra and intensity is greater than intensity for water z=2)
 - Orange/green shading indicates the extent of the peak, which is how far the algorithm goes out from the theoretical m/z peak to collect other spectra for the Gaussian curve orange if the peak is not significant, turns green if the fit is significant
 Plue shading is twice the extent (orange/green shading) which is the range of m/z that the
 - O Blue shading is twice the extent (orange/green shading), which is the range of m/z that the fitting algorithm looks in only visible if range is within the x-axis range
- Title includes the file name (without .mzML.gz) and the precursor loss composite type, then absolute value of the difference between the phosphoric acid z=2 and water loss z=2 delta m/z's in plots (1) and (2), and the ratio of phosphoric acid z=2 intensity to water loss z=2 intensity in the fitted peak (red) in plot (2)
 - Absolute difference between delta m/z's closer to 0 indicates lower error and if the fitted
 peaks are offset from the theoretical m/z (vertical line), this is likely due to absolute calibration
 offset
 - Intensity ratio is found by measuring the amplitude minus floor of the fitted Gaussian peak for phosphoric acid z=2, divided by that for water loss z=2